

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-BCL6 (Ser333) Polyclonal Antibody

多克隆抗体

产品货号: PA01710

第三版

描述:Transcriptional repressor which is required for germinal center formation and antibody affinity maturation. Probably plays an important role in lymphomagenesis.

抗原:Synthesized peptide derived from human Bcl-6 around the phosphorylation site of Ser333

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in germinal center T and B cells and in primary immature dendritic cells.