

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-SYT (Ser309/306) Polyclonal Antibody

多克隆抗体

产品货号: PA01756

第三版

描述: May have a regulatory role in the membrane interactions during trafficking of synaptic vesicles at the active zone of the synapse. It binds acidic phospholipids with a specificity that requires the presence of both an acidic head group and a diacyl backbone.

抗原: Synthesized peptide derived from human Synaptotagmin 1/2 around the phosphorylation site of Ser309/306

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: SYT1 is Expressed in melanocytes (PubMed:23999003). SYT2 is Expressed in melanocytes (PubMed:23999003).