

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Phospho-STK3/4 (Thr183) Polyclonal Antibody

多克隆抗体

产品货号: PA01803

第三版

描述:STK3 (Serine/Threonine Kinase 3) is a Protein Coding gene. Among its related pathways are Apoptosis and Autophagy and Sweet Taste Signaling. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is STK4. STK4 (Serine/Threonine Kinase 4) is a Protein Coding gene. Diseases associated with STK4 include T-Cell Immunodeficiency, Recurrent Infections, Autoimmunity, And Cardiac Malformations and Hyper-Ige Recurrent Infection Syndrome. Among its related pathways are Apoptosis and Autophagy and Glioma. GO annotations related to this gene include protein homodimerization activity and transferase activity, transferring phosphorus-containing groups. An important paralog of this gene is STK3.

抗原:Synthesized peptide derived from human Krs-1/2 around the phosphorylation site of Thr183

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: