

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-OPRM1 (Ser375) Polyclonal Antibody

多克隆抗体

产品货号: PA01877

第三版

描述: This gene encodes one of three opioid receptors. The mu opioid receptor is the principal target of endogenous opioid peptides and opioid analgesic agents such as beta-endorphin and enkephalins. The NM_001008503.1:c.118A>G allele had been associated with opioid and alcohol addiction and variations in pain sensitivity but evidence is conflicting. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Synthesized peptide derived from human MOR-1 around the phosphorylation site of Ser375

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed in brain. Isoform 16 and isoform 17 are detected in brain.