

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PAFAH Polyclonal Antibody

多克隆抗体

产品货号: PA01934

第三版

描述:PAFAH(platelet-activating factor acetylhydrolase, plasma) is also named as PLA2G7(phospholipase A2, group VII), Lp-PLA2 (Lipoprotein-associated phospholipase A2), LDL-PLA2, LDL-PLA and belongs to the AB hydrolase superfamily and lipase family. It is an enzyme which catalyzes the hydrolysis of acetyl ester at the sn-2 position of PAF. The family of PAF-Ahs consists of two intracellular isoforms (Ib and II), and one secreted isoform (plasma). Plasma PAF-AH is a monomeric polypeptide with a molecular weight of 45 kDa and it can show broad bands with a molecular weight of 43-67 kDa that it can be extensively N-glycosylated(PMID:12547653). Defects in PLA2G7 are the cause of platelet-activating factor acetylhydrolase deficiency (PAFAD), susceptibility to asthma (ASTHMA), susceptibility to atopic hypersensitivity (ATOPY).

抗原:Recombinant Human Platelet-activating factor acetylhydrolase protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: