

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse IL1A Polyclonal Antibody

多克隆抗体

产品货号: PA01985

第三版

描述:proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release.

抗原:Recombinant Human Interleukin-1 alpha protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: