

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Rat CD86 Polyclonal Antibody

多克隆抗体

产品货号: PA01992

第三版

**描述:** Receptor involved in the costimulatory signal essential for T-lymphocyte proliferation and interleukin-2 production, by binding CD28 or CTLA-4. May play a critical role in the early events of T-cell activation and costimulation of naive T-cells, such as deciding between immunity and anergy that is made by T-cells within 24 hours after activation. Isoform 2 interferes with the formation of CD86 clusters, and thus acts as a negative regulator of T-cell activation.

**抗原:** Recombinant human T-lymphocyte activation antigen CD86 protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Expressed by activated B-lymphocytes and monocytes.