

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ASCC1 Polyclonal Antibody

多克隆抗体

产品货号: PA01994

第三版

描述:Enhances NF-kappa-B, SRF and AP1 transactivation. In cells responding to gastrin-activated paracrine signals, it is involved in the induction of SERPINB2 expression by gastrin. May also play a role in the development of neuromuscular junction.

抗原:Recombinant Human Activating signal cointegrator 1 complex subunit 1 protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.