



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Rat IκB alpha Polyclonal Antibody

多克隆抗体

产品货号: PA01995

第三版

描述:Inhibits the activity of dimeric NF-kappa-B/REL complexes by trapping REL dimers in the cytoplasm through masking of their nuclear localization signals. On cellular stimulation by immune and proinflammatory responses, becomes phosphorylated promoting ubiquitination and degradation, enabling the dimeric RELA to translocate to the nucleus and activate transcription.

抗原:Recombinant human NF-kappa-B inhibitor alpha protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: