

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse TGFA Polyclonal Antibody

多克隆抗体

产品货号: PA02006

第三版

描述: TGF alpha is a mitogenic polypeptide that is able to bind to the EGF receptor/EGFR and to act synergistically with TGF beta to promote anchorage-independent cell proliferation in soft agar.

抗原: Recombinant human Protransforming growth factor alpha protein

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Isoform 1, isoform 3 and isoform 4 are expressed in keratinocytes and tumor-derived cell lines.