



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse PDGFRA Polyclonal Antibody

多克隆抗体

产品货号: PA02023

第三版

描述: Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

抗原: Recombinant human PDGFRA protein

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in primary and metastatic colon tumors and in normal colon tissue. Tumors may express a different isoform to that found in normal tissue.