
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ACER3 Polyclonal Antibody

多克隆抗体

产品货号: PA02053

第三版

描述:ACER3 (Alkaline Ceramidase 3) is a Protein Coding gene. Diseases associated with ACER3 include Farber Lipogranulomatosis and Partial Fetal Alcohol Syndrome. Among its related pathways are Sphingolipid metabolism and Metabolism. GO annotations related to this gene include hydrolase activity, acting on carbon-nitrogen (but not peptide) bonds, in linear amides and phytoceramidase activity.

抗原:Synthesized peptide derived from the Internal region of human APHC

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Ubiquitously expressed. Highest expression in placenta.