

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOF Polyclonal Antibody

多克隆抗体

产品货号: PA02066

第三版

描述:The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol. APOF (Apolipoprotein F) is a Protein Coding gene. Among its related pathways are Lipoprotein metabolism and Metabolism. GO annotations related to this gene include receptor binding and lipid transporter activity.

抗原:Synthesized peptide derived from the C-terminal region of human ApoF

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed by the liver and secreted in plasma.