

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AQP12 Polyclonal Antibody

多克隆抗体

产品货号: PA02072

第三版

**描述:**AQP12A (Aquaporin 12A) is a Protein Coding gene. Among its related pathways are Aquaporin-mediated transport and Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds. GO annotations related to this gene include transporter activity. An important paralog of this gene is AQP12B.AQP12B (Aquaporin 12B) is a Protein Coding gene. Among its related pathways are Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include transporter activity. An important paralog of this gene is AQP12A.

**抗原:**Synthesized peptide derived from the C-terminal region of human AQP12

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Restricted to the pancreas.