

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ARX Polyclonal Antibody

多克隆抗体

产品货号: PA02102

第三版

描述: This gene is a homeobox-containing gene expressed during development. The expressed protein contains two conserved domains, a C-peptide (or aristaless domain) and the prd-like class homeobox domain. It is a member of the group-II aristaless-related protein family whose members are expressed primarily in the central and/or peripheral nervous system. This gene is thought to be involved in CNS development. Expansion of a polyalanine tract and other mutations in this gene cause X-linked mental retardation and epilepsy. ARX (Aristaless Related Homeobox) is a Protein Coding gene. Diseases associated with ARX include Lissencephaly, X-Linked 2 and Epileptic Encephalopathy, Early Infantile, 1. Among its related pathways are Ectoderm Differentiation. GO annotations related to this gene include chromatin binding and transcriptional repressor activity, RNA polymerase II distal enhancer sequence-specific binding. An important paralog of this gene is RAX.

抗原: Synthesized peptide derived from the Internal region of human ARX

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed predominantly in fetal and adult brain and skeletal muscle. Expression is specific to the telencephalon and ventral thalamus. There is an absence of expression in the cerebellum throughout development and also in adult.