

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ADIPOQ Polyclonal Antibody

多克隆抗体

产品货号: PA02113

第三版

**描述:** Rat gAcrp30 (amp-1) is a 16.6 kDa protein that corresponds to the C-terminal globular domain of adiponectin (amp-1). gAcrp30/Adipolean is expressed and secreted exclusively by adipose tissue. This protein possesses unique pharmacological properties for regulation of body weight and lipid metabolism. It is a key component in controlling energy homeostasis. Recently, gAcrp30/Adipolean has been shown to inhibit the proliferation of myelomonocytic progenitors. The Recombinant rat gAcrp30/Adipolean is a 16.6 kDa protein consisting of 145 amino acid residues

**抗原:** Synthesized peptide derived from the N-terminal region of human Adiponectin

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Synthesized exclusively by adipocytes and secreted into plasma.