
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Cleaved-CASP2 p18 (G170) Polyclonal Antibody

多克隆抗体

产品货号: PA02118

第三版

描述:Involved in the activation cascade of caspases responsible for apoptosis execution. Might function by either activating some proteins required for cell death or inactivating proteins necessary for cell survival.

抗原:Synthesized peptide derived from the Internal region of human Caspase-2 p18

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed at higher levels in the embryonic lung, liver and kidney than in the heart and brain. In adults, higher level expression is seen in the placenta, lung, kidney, and pancreas than in the heart, brain, liver and skeletal muscle.