

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey Cleaved-CTSZ (L62) Polyclonal Antibody

多克隆抗体

产品货号: PA02133

第三版

描述:The protein encoded by this gene is a lysosomal cysteine proteinase and member of the peptidase C1 family. It exhibits both carboxy-monopeptidase and carboxy-dipeptidase activities. The encoded protein has also been known as cathepsin X and cathepsin P. This gene is expressed ubiquitously in cancer cell lines and primary tumors and, like other members of this family, may be involved in tumorigenesis.

抗原:Synthesized peptide derived from the Internal region of human Cathepsin Z

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed.