

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MAP3K5 Polyclonal Antibody

多克隆抗体

产品货号: PA02157

第三版

描述:Component of a protein kinase signal transduction cascade. Phosphorylates and activates MAP2K4 and MAP2K6, which in turn activate the JNK and p38 MAP kinases, respectively. Overexpression induces apoptotic cell death.

抗原:Synthesized peptide derived from human ASK 1 around the non-phosphorylation site of Ser83.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Abundantly expressed in heart and pancreas.