

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ATXN2L Polyclonal Antibody

多克隆抗体

产品货号: PA02164

第三版

描述: This gene encodes an ataxin type 2 related protein of unknown function. This protein is a member of the spinocerebellar ataxia (SCAs) family, which is associated with a complex group of neurodegenerative disorders. Several alternatively spliced transcripts encoding different isoforms have been found for this gene. ATXN2L (Ataxin 2 Like) is a Protein Coding gene. GO annotations related to this gene include poly(A) RNA binding. An important paralog of this gene is ATXN2.

抗原: Synthesized peptide derived from the Internal region of human Ataxin-2L

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed at high levels in thymus, lymph node, spleen, fetal kidney and adult testis. Constitutively associated with MPL and EPOR in hematopoietic cells.