

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BMP3 Polyclonal Antibody

多克隆抗体

产品货号: PA02246

第三版

**描述:**BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation. Negatively regulates bone density. Antagonizes the ability of certain osteogenic BMPs to induce osteoprogenitor differentiation and ossification. Expressed in adult and fetal cartilage. Highly expressed in fracture tissue, particularly in osteoblasts, osteoclasts and chondroblasts.

**抗原:**Synthesized peptide derived from the Internal region of human BMP-3.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in adult and fetal cartilage.