

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL20 Polyclonal Antibody

多克隆抗体

产品货号: PA02257

第三版

描述:The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis.

抗原:Synthetic peptide of human IL-20

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed at very low levels in skin, trachea, and other tissues.