

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HMG20B Polyclonal Antibody

多克隆抗体

产品货号: PA02258

第三版

**描述:** Required for correct progression through G2 phase of the cell cycle and entry into mitosis. Required for RCOR1/CoREST mediated repression of neuronal specific gene promoters.

**抗原:** Synthesized peptide derived from the N-terminal region of human BRAF35

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Ubiquitously expressed in adult tissues.