

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BRMS1 Polyclonal Antibody

多克隆抗体

产品货号: PA02272

第三版

描述: This gene reduces the metastatic potential, but not the tumorigenicity, of human breast cancer and melanoma cell lines. The protein encoded by this gene localizes primarily to the nucleus and is a component of the mSin3a family of histone deacetylase complexes (HDAC). The protein contains two coiled-coil motifs and several imperfect leucine zipper motifs. Alternative splicing results in two transcript variants encoding different isoforms. BRMS1 (Breast Cancer Metastasis Suppressor 1) is a Protein Coding gene. Diseases associated with BRMS1 include Breast Cancer. Among its related pathways are Activated PKN1 stimulates transcription of AR (androgen receptor) regulated genes KLK2 and KLK3 and Chromatin organization. GO annotations related to this gene include histone deacetylase binding and histone deacetylase activity. An important paralog of this gene is BRMS1L.

抗原: Synthesized peptide derived from the Internal region of human BRMS-1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expression levels are higher in term placentas than in early placentas. Low levels of expression observed in normal pregnancies and in molar pregnancies.