

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BRS3 Polyclonal Antibody

多克隆抗体

产品货号: PA02276

第三版

**描述:**The protein encoded by this gene is a G protein-coupled membrane receptor that binds bombesin-like peptides. This binding results in activation of a phosphatidylinositol-calcium second messenger system, with physiological effects including regulation of metabolic rate, glucose metabolism, and hypertension.

**抗原:**Synthesized peptide derived from the Internal region of human BRS-3

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**In germ cells in testis. Lung carcinoma cells.