

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IGFBP6 Polyclonal Antibody

多克隆抗体

产品货号: PA02324

第三版

描述:IGFBP6 (Insulin Like Growth Factor Binding Protein 6) is a Protein Coding gene. Diseases associated with IGFBP6 include Exfoliative Ichthyosis and Insulin-Like Growth Factor I. Among its related pathways are Metabolism of proteins and Development IGF-1 receptor signaling. GO annotations related to this gene include insulin-like growth factor binding and insulin-like growth factor II binding. An important paralog of this gene is IGFBP5.

抗原:Synthesized peptide derived from the Internal region of human IGFBP6.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: