

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human XK Polyclonal Antibody

多克隆抗体

产品货号: PA02331

第三版

描述: This locus controls the synthesis of the Kell blood group 'precursor substance' (Kx). Mutations in this gene have been associated with McLeod syndrome, an X-linked, recessive disorder characterized by abnormalities in the neuromuscular and hematopoietic systems. The encoded protein has structural characteristics of prokaryotic and eukaryotic membrane transport proteins.

抗原: Synthetic peptide of human XK

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: High levels in skeletal muscle, heart, brain, and pancreas; low levels in placenta, lung, liver, and kidney.