

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human XKRX Polyclonal Antibody

多克隆抗体

产品货号: PA02344

第三版

描述: This gene encodes a protein that is related to a component of the XK/Kell complex of the Kell blood group system. The encoded protein includes several transmembrane domains, is known to be exposed to the cell surface, and may function as a membrane transporter.

抗原: Synthetic peptide of human XKRX

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed predominantly in the placenta, in syncytiotrophoblasts. Moderate levels in the adrenal gland, low levels in the trachea and very low levels in the bone marrow.