

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GORAB Polyclonal Antibody

多克隆抗体

产品货号: PA02345

第三版

描述: This gene encodes a member of the golgin family, a group of coiled-coil proteins localized to the Golgi. The encoded protein may function in the secretory pathway. The encoded protein, which also localizes to the cytoplasm, was identified by interactions with the N-terminal kinase-like protein, and thus it may function in mitosis. Mutations in this gene have been associated with geroderma osteodysplastica. Alternatively spliced transcript variants have been described. GORAB (Golgin, RAB6 Interacting) is a Protein Coding gene. Diseases associated with GORAB include Geroderma Osteodysplasticum and Geroderma Osteodysplastica. Among its related pathways are DNA Damage Response.

抗原: Synthesized peptide derived from the N-terminal region of human SCYL1BP1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: