

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZNF473 Polyclonal Antibody

多克隆抗体

产品货号: PA02364

第三版

描述:Involved in histone 3'-end pre-mRNA processing by associating with U7 snRNP and interacting with SLBP/pre-mRNA complex. Increases histone 3'-end pre-mRNA processing but has no effect on U7 snRNP levels, when overexpressed. Required for cell cycle progression from G1 to S phases.

抗原:Synthetic peptide of human ZNF473

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: