
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CA13 Polyclonal Antibody

多克隆抗体

产品货号: PA02365

第三版

描述:Carbonic anhydrases (CAs) are a family of zinc metalloenzymes. For background information on the CA family, see MIM 114800.CA13 (Carbonic Anhydrase 13) is a Protein Coding gene. Among its related pathways are Nitrogen metabolism and Metabolism. GO annotations related to this gene include carbonate dehydratase activity. An important paralog of this gene is CA2.

抗原:Synthesized peptide derived from the Internal region of human CA XIII

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in thymus, small intestine, spleen, prostate, ovary, colon and testis.