

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat COQ8A Polyclonal Antibody

多克隆抗体

产品货号: PA02368

第三版

**描述:** This gene encodes a mitochondrial protein similar to yeast ABC1, which functions in an electron-transferring membrane protein complex in the respiratory chain. It is not related to the family of ABC transporter proteins. Expression of this gene is induced by the tumor suppressor p53 and in response to DNA damage, and inhibiting its expression partially suppresses p53-induced apoptosis. Alternatively spliced transcript variants have been found; however, their full-length nature has not been determined. COQ8A (Coenzyme Q8A) is a Protein Coding gene. Diseases associated with COQ8A include Coenzyme Q10 Deficiency, Primary, 4 and Coq8a-Related Coenzyme Q10 Deficiency. An important paralog of this gene is COQ8B.

**抗原:** Synthesized peptide derived from the Internal region of human CABC1.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Ubiquitously expressed with a relatively greater abundance in heart and skeletal muscle.