

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CABP7 Polyclonal Antibody

多克隆抗体

产品货号: PA02379

第三版

**描述:**CABP7 (Calcium Binding Protein 7) is a Protein Coding gene. Among its related pathways are Visual Cycle in Retinal Rods. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CALN1.

**抗原:**Synthesized peptide derived from the N-terminal region of human CaBP7

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**