

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Cadherin-19 Polyclonal Antibody

多克隆抗体

产品货号: PA02389

第三版

**描述:** This gene is one of three related type II cadherin genes situated in a cluster on chromosome 18. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein containing five extracellular cadherin repeats. Loss of cadherins may be associated with cancer formation. Alternative splicing results in multiple transcript variants for this gene.

**抗原:** Synthesized peptide derived from the Internal region of human Cadherin-19

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Expressed in many tissues, with the exception of uterus.