

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Cadherin-like protein 26 Polyclonal Antibody

多克隆抗体

产品货号: PA02390

第三版

描述: This gene encodes a member of the cadherin protein family. Cadherins are a family of calcium-dependent adhesion molecules that mediate cell-cell adhesion in all solid tissues and modulate a wide variety of processes, including cell polarization, migration and differentiation. Cadherin domains occur as repeats in the extracellular region and are thought to contribute to the sorting of heterogeneous cell types and the maintenance of orderly structures such as epithelium. This protein is expressed in gastrointestinal epithelial cells and may be upregulated during allergic inflammation. This protein interacts with alpha integrins and may also be involved in leukocyte migration and adhesion. CDH26 (Cadherin 26) is a Protein Coding gene. Among its related pathways are Nanog in Mammalian ESC Pluripotency and ERK Signaling. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CDH2.

抗原: Synthesized peptide derived from the N-terminal region of human Cadherin-26

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: