

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CAPN12 Polyclonal Antibody

多克隆抗体

产品货号: PA02404

第三版

描述:The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes a member of the calpain large subunit family. CAPN12 (Calpain 12) is a Protein Coding gene. Diseases associated with CAPN12 include Autosomal Recessive Congenital Ichthyosis. Among its related pathways are Degradation of the extracellular matrix and Apoptosis Pathway. GO annotations related to this gene include calcium ion binding and calcium-dependent cysteine-type endopeptidase activity. An important paralog of this gene is CAPN1.

抗原:Synthesized peptide derived from the Internal region of human Calpain 12

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: