
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CASP10 B/C Polyclonal Antibody

多克隆抗体

产品货号: PA02452

第三版

描述:Involved in the activation cascade of caspases responsible for apoptosis execution. Recruited to both Fas- and TNFR-1 receptors in a FADD dependent manner. May participate in the granzyme B apoptotic pathways. Cleaves and activates caspase-3, -4, -6, -7, -8, and -9. Hydrolyzes the small- molecule substrates, Tyr-Val-Ala-Asp-AMC and Asp-Glu-Val-Asp-AMC. Isoform C is proteolytically inactive.

抗原:Synthesized peptide derived from the C-terminal region of human Caspase-10 B/C

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detectable in most tissues. Lowest expression is seen in brain, kidney, prostate, testis and colon.