

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TYMP Polyclonal Antibody

多克隆抗体

产品货号: PA02465

第三版

描述: May have a role in maintaining the integrity of the blood vessels. Has growth promoting activity on endothelial cells, angiogenic activity in vivo and chemotactic activity on endothelial cells in vitro. Catalyzes the reversible phosphorolysis of thymidine. The produced molecules are then utilized as carbon and energy sources or in the rescue of pyrimidine bases for nucleotide synthesis.

抗原: Synthesized peptide derived from the N-terminal region of human PD-ECGF.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: