
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR87/95 Polyclonal Antibody

多克隆抗体

产品货号: PA02486

第三版

描述: Receptor for lysophosphatidic acid (LPA). Necessary for p53/TP53-dependent survival in response to DNA damage.

抗原: Synthesized peptide derived from the N-terminal region of human GPR87/95.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in placenta and prostate. Weaker expression in thymus. Not expressed in thalamus, hippocampus, pons or cerebellum. Overexpressed in squamous cell carcinoma of the lung.