

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CTSG Polyclonal Antibody

多克隆抗体

产品货号: PA02489

第三版

描述: Serine protease with trypsin- and chymotrypsin-like specificity. Cleaves complement C3. Has antibacterial activity against the Gram-negative bacterium *P.aeruginosa*, antibacterial activity is inhibited by LPS from *P.aeruginosa*, Z-Gly-Leu-Phe-CH₂Cl and phenylmethylsulfonyl fluoride.

抗原: Synthesized peptide derived from the Internal region of human Cathepsin G.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: