

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NCBP2 Polyclonal Antibody

多克隆抗体

产品货号: PA02513

第三版

描述:NCBP2 is a component of the nuclear cap-binding protein complex (CBC), which binds to the monomethylated 5' cap of nascent pre-mRNA in the nucleoplasm. It has an RNP domain commonly found in RNA binding proteins, and contains the cap-binding activity.

抗原:Synthesized peptide derived from the N-terminal region of human CBP20

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: