

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CD130 Polyclonal Antibody

多克隆抗体

产品货号: PA02530

第三版

**描述:**Signal-transducing molecule. The receptor systems for IL6, LIF, OSM, CNTF, IL11, CTF1 and BSF3 can utilize gp130 for initiating signal transmission. Binds to IL6/IL6R (alpha chain) complex, resulting in the formation of high-affinity IL6 binding sites, and transduces the signal. Does not bind IL6. May have a role in embryonic development (By similarity). The type I OSM receptor is capable of transducing OSM-specific signaling events.

**抗原:**Synthesized peptide derived from human CD130 around the non-phosphorylation site of Ser782.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Found in all the tissues and cell lines examined. Expression not restricted to IL6 responsive cells.