

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LILRA1 Polyclonal Antibody

多克隆抗体

产品货号: PA02558

第三版

**描述:** This gene encodes an activating member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein is predominantly expressed in B cells, interacts with major histocompatibility complex class I ligands, and contributes to the regulation of immune responses. Alternative splicing results in multiple transcript variants encoding different isoforms.

**抗原:** Synthesized peptide derived from the N-terminal region of human LIR-6.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Detected in monocytes and B-cells.