

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CD110 Polyclonal Antibody

多克隆抗体

产品货号: PA02562

第三版

描述:The protein encoded by the c-mpl gene, CD110, is a 635 amino acid transmembrane domain, with two extracellular cytokine receptor domains and two intracellular cytokine receptor box motifs. TPO-R deficient mice were severely thrombocytopenic, emphasizing the important role of CD110 and thrombopoietin in megakaryocyte and platelet formation. Upon binding of thrombopoietin CD110 is dimerized and the JAK family of non-receptor tyrosine kinases, as well as the STAT family, the MAPK family, the adaptor protein Shc and the receptors themselves become tyrosine phosphorylated. MPL (MPL Proto-Oncogene, Thrombopoietin Receptor) is a Protein Coding gene. Diseases associated with MPL include Thrombocytopenia, Congenital Amegakaryocytic and Myelofibrosis With Myeloid Metaplasia, Somatic. Among its related pathways are Response to elevated platelet cytosolic Ca²⁺ and NF-kappaB Signaling. GO annotations related to this gene include transmembrane signaling receptor activity and thrombopoietin receptor activity. An important paralog of this gene is EPOR.

抗原:Synthesized peptide derived from the Internal region of human CD110.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Expressed at a low level in a large number of cells of hematopoietic origin. Isoform 1 and isoform 2 are always found to be coexpressed.