



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD3EAP Polyclonal Antibody

多克隆抗体

产品货号: PA02570

第三版

描述:CD3EAP (CD3e Molecule Associated Protein) is a Protein Coding gene. Among its related pathways are RNA Polymerase I Promoter Escape and Activated PKN1 stimulates transcription of AR (androgen receptor) regulated genes KLK2 and KLK3. GO annotations related to this gene include poly(A) RNA binding and DNA-directed 5-3 RNA polymerase activity.

抗原:Synthesized peptide derived from the C-terminal region of human CD3EAP

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: