

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD302 Polyclonal Antibody

多克隆抗体

产品货号: PA02575

第三版

**描述:**CD302 is a C-type lectin receptor involved in cell adhesion and migration, as well as endocytosis and phagocytosis. CD302 (CD302 Molecule) is a Protein Coding gene. Diseases associated with CD302 include Hodgkin Lymphoma. GO annotations related to this gene include carbohydrate binding.

**抗原:**Synthesized peptide derived from the Internal region of human DCL-1.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in myeloid and B lymphoid cell lines. Isoform 2 and isoform 3 are expressed in malignant Hodgkin lymphoma cells called Hodgkin and Reed-Sternberg (HRS) cells.