

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CMTM5 Polyclonal Antibody

多克隆抗体

产品货号: PA02596

第三版

描述: This gene belongs to the chemokine-like factor gene superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. Highly expressed in the brain.

抗原: Synthesized peptide derived from the N-terminal region of human CMTM5.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in the brain.