

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse POLR3G Polyclonal Antibody

多克隆抗体

产品货号: PA02600

第三版

描述:POLR3G (RNA Polymerase III Subunit G) is a Protein Coding gene. Diseases associated with POLR3G include Leukemia, Acute Myelomonocytic, Somatic, Somatic. Among its related pathways are Innate Immune System and RNA Polymerase III Transcription Initiation. GO annotations related to this gene include DNA-directed 5-3 RNA polymerase activity and RNA polymerase III activity. An important paralog of this gene is POLR3GL.

抗原:Synthesized peptide derived from the Internal region of human POLR3G.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Barely detectable in differentiated tissues. Expressed in embryonic stem cells and in other dividing cells, such as some tumor cell lines.