

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SAMHD1 Polyclonal Antibody

多克隆抗体

产品货号: PA02603

第三版

描述: This gene may play a role in regulation of the innate immune response. The encoded protein is upregulated in response to viral infection and may be involved in mediation of tumor necrosis factor- α proinflammatory responses. Mutations in this gene have been associated with Aicardi-Goutieres syndrome.

抗原: Synthesized peptide derived from the Internal region of human SAMHD1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in heart, skeletal muscle, spleen, liver, small intestine, placenta, lung and peripheral blood leukocytes. No expression is seen in brain and thymus.