

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat CDCA3 Polyclonal Antibody

多克隆抗体

产品货号: PA02633

---

第三版

**描述:**CDCA3 (Cell Division Cycle Associated 3) is a Protein Coding gene. Among its related pathways are Cell cycle\_Role of APC in cell cycle regulation.

**抗原:**Synthesized peptide derived from the C-terminal region of human CdcA3

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**